



User Manual

retinaWISE

General Information	
Filename:	02-1_UserManual.docx
Number of Pages:	38
Version number:	2.0
Issue Date:	Tuesday, March 05, 2024



Table of Content

1	GENERAL INFORMATION	5
1.1	Retina Wise overview	5
1.2	Manufacturer	5
2	SAFETY	6
2.1	General Safety Information	6
2.2	Training of the Staff	6
2.3	Working Area	6
2.4	Interference with Other Devices	6
2.5	Instructions for the Device Disposal	7
2.6	Labels	8
3	ENVIRONMENTAL CONDITIONS AND FIRST START-UP	10
3.1	Responsibility of the Customer	10
3.2	Requirements for the Electrical Connection	10
3.3	Temperature, Pressure and Humidity	10
3.4	Packaging	10
3.5	Physical Installation	10
3.6	Labelling Check	10
3.7	First start-up	10
3.8	Device Connections	12
3.9	Software installation	12
4	DESCRIPTION OF THE DEVICE	13
4.1	System Overview	13
4.2	Hardware	13
4.3	User Interface	14
5	OPERATING INSTRUCTIONS	17
5.1	First configuration	17
5.2	Connect the communication port	17
5.3	Connect the power supply and power on the Device	18
5.4	Start the software	18
5.5	Alignment procedure	19
5.6	Loading a protocol	23
5.7	Launch the experiment	23
5.8	Analyse the results	24
6	CLINICAL APPLICATION	26
6.1	Training requirements	26
6.2	Intended use	26

6.3 Patient Population	26
6.4 Contraindications	26
6.5 Side Effects	26
6.6 Residual Risk.....	26
6.7 Warnings related to the applications.....	26
6.8 Precautions.....	27
6.9 Adverse events management	27
7 MAINTENANCE	28
7.1 Checking the line cable	28
7.2 Checking the labels	28
7.3 Cleaning.....	28
8 TROUBLESHOOTING	30
9 CUSTOMER SERVICE	31
9.1 Warranty and Manufacturer's responsibilities	31
9.2 Repairs and modifications of the Device	31
9.3 Contacts	31
9.4 Reporting of serious incidents	32
10 TECHNICAL SPECIFICATIONS	33
10.1 General Specifications.....	33
10.2 LED specifications.....	33
11 GLOSSARY	34
12 APPENDIX A: EMC TABLES	35
13 APPENDIX B:	38

1 GENERAL INFORMATION

1.1 Retina Wise overview

The retinaWISE Medical Device evaluates the sensitivity of retina tissue to specific light stimuli. The evaluation allows the operators to gather objective information that will let them conclude on sensitivity and health of retina tissues. retinaWISE is a medical device for professional use only.



Warning: Any use of this Device deviating from the prescriptions included in this Manual is strictly forbidden and would be out of the liability of the Manufacturer

The functioning of the retinaWISE Device is based on the monitoring of patients' pupillary reflex by measuring the pupil's diameter dynamic as a function of scotopic light stimuli with defined spectroscopic characteristics.

1.2 Manufacturer

The retinaWise medical Device is manufactured following the REGULATION (EU) 2017/745.

Model		retinaWISE
Category of the Device		Medical Device with relative measuring function
Classification		Class I (According to REGULATION (EU) 2017/745)
	Manufactured by	Oculox Technologies SA Via Industria 3 CH-6933 Muzzano Switzerland Tel: +41 91 210 89 61
EC REP		Qarad EC-REP BV Pas 257 2440 Geel Belgium SRN: BE-AR-000000040

2 SAFETY

2.1 General Safety Information

0. To safely use the Device, knowing all the safety rules according to International Standards is necessary.
1. All the people operating this equipment must understand the operation and safety instructions specified in this Manual.
2. Only authorised individuals with appropriate medical training and knowledge should operate the system.
3. Only authorised Service personnel should access the internal/electrical components of the system.
4. The User Manual must be available in the working area of the Device.
5. Warning labels must be all kept in good condition.



CAUTION: *Using controls, adjustments, or procedures other than those specified herein may result in hazardous radiation exposure.*



CAUTION: *US Federal law restricts this Device to sale by, or on the order of, a physician/surgeon.*

2.2 Training of the Staff

retinaWISE shall be operated only by medical staff with experience in Ophthalmology. According to their expertise, such professionals can decide the correct use of the Device depending on the type of application (Optometrist, Ophthalmologist).

2.3 Working Area

The Device shall be placed on a flat and stable surface.

The Device shall be placed in a darkened room to allow scotopic measurements.

2.4 Interference with Other Devices

This Device is intended to be used in connection with a computer.

Minimum Computer Configuration/Requirements	
Processor (CPU)	Intel Core i3 or equivalent
Operating System	Microsoft Windows 10 x64
RAM Memory [GB]	8 GB RAM
Storage [GB]	100 GB internal storage drive
USB adapter	3.0 type A (or newer)
Monitor/Display	1920x1080 resolution
Antivirus – Software	Update

This Device can be disturbed by the interference of external electromagnetic fields generated by other electrical devices in the closest proximity of the Device.

The Device can work only with the approved software installed on a PC duly protected against malware through state-of-the-art protections such as antivirus, firewall etc.



CAUTION: *Mobile phones and similar electrical devices must be switched off when the Device is working.*



CAUTION: USB port is accessible only by trained and authorized service personnel. Restriction of use are present

This Device must be installed and used according to EMC (Electromagnetic Compatibility) information described in the tables reported in Appendix A: EMC Tables.

2.5 Instructions for the Device Disposal

At the end of its lifetime, the Device must be handled according to the National or Local regulations for electrical and electronic equipment disposal.

The Device is subject to national standards regulating waste disposal, such as electrical equipment.



It is forbidden to dispose of the Device as municipal waste. Instead, it must be collected separately according to the WEEE Directive (Waste Electrical and Electronic Equipment).

The penalties for violating law requirements are severe.

2.6 Labels

Please find a representative image of the official label applied to the Device. The labelling plan and the final solution are discussed exhaustively in the corresponding documents ([02-2 LabelingPlan](#)).

2.6.1 Main Device's label




REF

A001-007-v01


SN

ReW2401001





2024-01-19

UDI

(01) 7649988168413
(11) 240119
(21) ReW2401001



CE

100-240 VAC 50/60 Hz
24 VDC 40 VA

 Use Only with power supply Meanwell GSM90B24-P1M

**Oculox Technologies SA**
Via Industria 3
CH-6933 Muzzano

EC

REP

Qarad EC-REP BV
Pas 257
B-2440 Geel
SRN:BE-AR-0000000040

2.6.2 Label positioning



3 ENVIRONMENTAL CONDITIONS AND FIRST START-UP

3.1 Responsibility of the Customer

The user is responsible for placing the Device on a stable table in a clean environment (protected from water and dust).

3.2 Requirements for the Electrical Connection

The Device shall be connected to a power socket in compliance with the local electrical safety regulations. Following the Standard IEC/EN 60884-1, the shape and type of the cable's connector shall be adapted to the local power socket of the country where the Device will be installed.



Warning: The Device shall be installed and used in compliance with your country's national or local requirements.

3.3 Temperature, Pressure and Humidity

The Device requires a dry, low-dust area with adequate ventilation. Air conditioning is preferable but not essential. Refer to Section 10 of this manual for temperature, pressure and humidity ranges.

3.4 Packaging

The Device is shipped in a dedicated carton box. Packaging includes conformal protective foam. The Client is responsible for reviewing delivered goods, prepare and position the Device in the designed working area ready for installation.

3.5 Physical Installation

The device is constructed to be well balanced, and it is designed to be positioned on a standard ophthalmic table. retinaWISE is provided with the possibility to be rigidly fixed on the ophthalmic table.



CAUTION: retinaWise need to be stably positioned on an ophthalmic table. Consider to use the provided screws to fix the device on the ophthalmic table

3.6 Labelling Check

The user's responsibility is to verify the integrity and readability of the security labels placed on the Device. Damaged labels shall be replaced immediately in compliance with the Labelling shown in [§2.6](#).

3.7 First start-up

The first start-up of the Device must be performed after having read and understood the User Manual. Please do not hesitate to contact our service personnel see §9.3 in case of questions or doubts.



Warning: *Do not start using the Device without reading this Manual. The warranty does not cover damage that occurred before installation.*

The received material shall be inspected immediately upon arrival on the following terms:

- Administrative check:
Number of packages
Sizes and weights
- Technical check:
Packaging condition

These checks must be made visually, with the most outstanding possible care and in the presence of the delivery employee.

The installation procedure shall be performed each time the Device is installed for the first time or after being transported through cars, elevators, trucks, air cargo, etc.

During installation, the Device shall be checked for proper operation, and eventual malfunctions due to transportation must be corrected.

The installation procedure also includes a training course from the Distributor to the user concerning the use of the medical Device.

The first turn-on procedure may take several hours. During this time, access to the installation site is forbidden.

The packed materials shall be checked immediately upon their arrival, if possible, in the presence of the shipper's delivery employee, as follows:

- Open the packaging and put the Device in a proper site for a general check
- Perform the following operations for the general check:
 - No obvious damage occurred during the shipping
- Perform further controls or additional tests
- Perform training for the End User on the operation of the Device

3.8 Device Connections

The Device is delivered with its specific power supply; no other power supplies are allowed.

The Device's connection with the pc is made with the delivered USB cable or another USB cable with at least a 3.0 communication standard.

3.9 Software installation

The installation of the software and the setting of the parameter for the usage of the Device has to be done by a technical specialist from Oculox Technologies.

4 DESCRIPTION OF THE DEVICE

retinaWise general description and details of the main components and accessories.

4.1 System Overview

The retinaWise presents as a standard ophthalmic appliance mounted on a joystick and interfaces the patients through a chin rest. The joystick allows the Operator to align the Device's projection optics to the patient's eyes. The Operator drives the measurements from the computer adequately connected to the Device.



4.2 Hardware

The Device is composed of an ophthalmic joystick that moves the optical system to allow the operators to align it with the patient eyes.

The Device is connected to a PC through a USB cable.

A power switch ON and OFF is mounted on the Device.

The software needs to be installed on the PC to allow the connection to the Device.

4.3 User Interface


The software that drives the Device is property of Oculox Technologies



CAUTION: Only the supplied software can be used to control the Device.

The user interface is composed of the following main zones:

1. Menu Bar
2. Patient's eye image stream
3. Summary of the measurement session
4. Graphic representation of excitation & data plots
5. Functional Button Zone



4.3.1 Menu Bar

In the menu bar, there are six extendable menus:

- Action
 - Users can start all the steps needed to perform the analysis; these commands are also accessible from the *Function Button Zone* in the main window.
- Edit
 - Users to edit the Device's settings; the section is reserved only for qualified technical personnel; in case of troubleshooting, contact Oculox Technologies
- Statistics
 - Device usage and pupil diameter statistics.
- Performance
 - Acquisition and data storage performances and details
- Tracking Status
 - Details and objective evaluations on the pupil detection algorithm
- Help
 - Software's version
 - Technical support - Automatic ticketing system in direct connection with technical support.


4.3.2 Patient's eye image stream

Realtime images are acquired from the monitoring cameras embedded into the Device; the user visually checks the alignment status to apply corrections by acting on the joystick to allow proper data collection.

To support the alignment procedure, an augmented reality algorithm projects the detected pupil circumference on top of the patient's eye streams.

4.3.3 Summary of a measurement session

Simple and synthetic portion providing a summary of patient data and the protocol's details.



4.3.4 Graphic representation of excitation & data plots

Excitation settings and data plots are presented in two distinct zones:

1. The two graphs on the higher portion plot the time information about the light stimuli
 - a. Left eye, and
 - b. Right eye
2. In the lower portion, the graph plots the real-time measure pupil's diameter for both eyes.

NOTE: the user can perform the same protocol several times; in this case, the previous results will be displayed in transparency on the lower graph.

4.3.5 Functional Button Zone

Control zone to successfully perform the analysis.

From left to right, the user will need to pass through every button which corresponds to different action:

- Aligning the patient
- Load the protocol
- Start the sequence of analysis
- Stop the sequence of analysis
 - Only if necessary

5 OPERATING INSTRUCTIONS


5.1 First configuration

Before starting to use the Device, please ensure that the Device is adequately placed on a flat, stable surface.

Oculox Technicians will do the initial configuration; the user does not need to take care of setting up the Device.


5.2 Connect the communication port

Connect the supplied USB cable via the supplied optoisolator, to the Device on one side and the computer on the other.





5.3 Connect the power supply and power on the Device

Connect the supplied power supply to the Device and position the main ON-OFF switch.



5.4 Start the software

Start the installed software retinaWISE "Run RetinaWISE or the executable file rewise.exe.




The programs start, and the Device is ready to operate.


5.5 Alignment procedure

The Device is built to be adjusted along multiple mechanical degrees of freedom. The strategy allows the alignment of a fixed patient on the moving Device's optics.

The joystick allows the horizontal movement to align with the right eye; by rotating the handle, the vertical position of the Device will be adapted to the patient's height.



In addition to the degrees of freedom provided by the joystick, it reveals to be vital the possibility to adapt the height of the headrest to adjust to the different head anatomies



Before starting the procedure, place the vertical position of the joystick in the middle position by rotating the handle.

5.5.1 Pre-alignment of the patient

Before inviting the patient to accommodate on Device's chin rest and avoid uncomfortable neck and back positions for the patient, adapt sequentially:

- i) The ophthalmic table height, and
- ii) The patient's chair height.

Once the patient is comfortably placed, invite him to put his chin on the chin guard and rest the forehead against the support.


Verify if the eyes are at the same level comparing with the camera of the Device, if not turn the handle on the headrest to adjust the vertical position of the chin guard.




5.5.2 Patient right eye

On the UI press the "Align Patient" button to start the alignment procedure.

Move the joystick horizontally and vertically (moving and rotating the handle) aligning the right eye of the patient in front of the Device by looking at the streamed image on the User Interface Window.




Move the Device to reach the situation where the small yellow circle is centred with the right pupil.




5.5.3 Patient left eye

Rotate the nut placed on the right side of the Device to adjust the ocular distance until the left eyes is detected by the software.



NOTE: It will be possible that the medics need to ask the patient to tilt a bit the head to align the right eye and the left eye horizontally.


The final alignment situation will be:



5.6 Loading a protocol


We assume that a protocol already exists. It can be loaded either with the button on the bottom or while the menu action.

Select file of the desired protocols from the opened window



5.7 Launch the experiment

After selected the protocol to use, press "Start Sequence" button. The Experiment is now performing, the window will be filled with the results of the experiment in real time.




5.8 Analyse the results

5.8.1 Application information

The results of the experiments are reported directly on the window, the two upper graphs show the light protocol directed on the eyes as function of the time for the left and right eye.

The lower graph shows the pupil diameter reaction to these light stimuli as a function of time.



5.8.2 Further analysis – saved data layout

The software saves automatically some files related to the experiment done under a predefined folders that can be chosen in according with the specialist that install the Device.

The data are exported to a .csv file with the following definition:

Time	Excitation index	Left - Is four	Left - Size M	Left - Area M	Left - Radius I	Left - Radius B	Left - Pos X	Left - Pos Y	Left - Distance from foci	Left - Leads ref	Right - Is four	Right - Size M	Right - Area M	Right - Radius I	Right - Radius B	Right - Pos X	Right - Pos Y	Right - Distance from foci	Right - Leads ref
2.002	2	VRAI	6.49392	33.1203	123.03398	168.2545	521.84285	553.56403	22.66434	24.7392	VRAI	6.41042	32.27479	128.55362	156.90578	475.4793	551.9385	64.81723	-0.50195
4.012	3	VRAI	6.51685	33.3553	124.75602	168.23983	521.51727	553.85495	23.08891	24.70817	VRAI	6.41012	32.30788	128.5278	156.71103	475.7285	552.08944	65.55229	-0.50195
6.002	4	VRAI	6.52891	33.47475	124.68025	167.62285	520.72183	553.55964	23.56525	24.71595	VRAI	6.43395	32.52011	128.57574	156.76987	475.8386	552.23778	64.56211	-0.50195
8.002	5	VRAI	6.50994	33.28464	123.81031	168.2342	522.03174	553.55511	22.50777	24.70428	VRAI	6.40571	32.26764	128.0737	156.2462	474.85245	551.64813	65.3826	-0.50195
10	5	VRAI	6.50951	33.28599	123.8441	168.17552	521.72455	553.83325	22.98739	24.70817	VRAI	6.41845	32.55563	123.05364	156.69783	475.71809	552.16774	65.0321	-0.50195
12	5	VRAI	6.49622	33.14453	124.4276	167.84205	521.50323	553.76111	23.05425	24.71006	VRAI	6.4166	32.33701	128.69637	157.04321	475.78934	552.06258	64.5503	-0.50195
14.02	6	VRAI	6.50246	33.2082	123.48027	168.11177	521.91705	554.26074	23.02359	24.71895	VRAI	6.40096	32.17859	128.45462	156.57082	475.4955	551.79547	64.8552	-0.50195
16	7	VRAI	6.49013	33.08241	123.57768	166.9056	521.41053	553.04359	22.71719	24.71036	VRAI	6.43753	32.54896	123.66443	156.80545	475.46668	551.9306	64.82302	-0.50195
18	7	VRAI	6.5052	33.23615	123.55039	168.1636	521.52155	554.02502	23.1896	24.71804	VRAI	6.41955	32.28711	128.59418	156.35564	475.3602	552.057	64.96438	-0.50195
20.01	8	VRAI	6.50342	33.21796	123.55473	168.03281	521.66431	553.32779	23.02529	24.71895	VRAI	6.41254	32.2361	128.86551	156.63666	475.6755	552.86327	64.73591	-0.50195
22	9	VRAI	6.50936	33.23879	123.71862	167.8933	521.33238	553.68445	23.16228	24.71894	VRAI	6.41822	32.32194	128.58274	157.10034	475.6239	552.08104	64.71322	-0.50195
24	9	VRAI	6.51853	33.37253	124.14056	168.01785	521.33038	553.78845	23.1839	24.72763	VRAI	6.43729	32.54593	128.87665	156.67442	476.13431	552.46894	64.29613	-0.50195
26	10	VRAI	6.49536	33.07451	123.59867	168.05198	522.23415	553.80127	22.44329	24.70817	VRAI	6.41119	32.28251	128.53091	156.90509	475.52422	551.68542	64.38714	-0.50195
28.01	11	VRAI	6.5094	33.18734	123.88386	168.05445	521.083	554.14311	22.38389	24.71894	VRAI	6.4181	32.33382	128.66652	157.05239	475.59778	551.13176	64.70875	-0.50195
30	11	VRAI	6.52181	33.40777	124.07828	168.37033	521.1347	554.06362	23.53422	24.70817	VRAI	6.44343	32.60801	128.55846	157.30356	476.254	552.04102	64.08689	-0.50195
32.01	11	VRAI	6.50788	33.28366	123.44262	168.14621	521.88878	554.3538	23.11177	24.71892	VRAI	6.41402	32.31005	128.75864	156.73839	475.55335	552.03375	64.72689	-0.50195
34.01	11	VRAI	6.53507	33.54207	124.46802	168.5708	521.83382	554.20819	23.53385	24.71036	VRAI	6.43457	32.41734	123.3976	156.63339	476.73076	551.8819	63.34033	-0.50195
36	12	VRAI	6.51644	33.33273	123.51718	168.67264	521.59122	554.21521	23.2822	24.70428	VRAI	6.418	32.35114	123.4086	156.24512	475.73068	552.03601	64.59333	-0.50195
38	13	VRAI	6.50539	33.22818	124.07182	167.42391	521.55646	553.52405	22.8706	24.72374	VRAI	6.42474	32.41915	123.59805	156.30156	475.71814	552.10144	64.63044	-0.50195
40	14	VRAI	6.50553	33.2313	123.71417	167.82189	521.71021	553.93072	23.02731	24.72763	VRAI	6.41553	32.32623	123.11497	156.47984	475.33021	551.45335	64.2351	-0.50195
42	15	VRAI	6.50584	33.24067	123.71534	167.58239	521.57825	554.0083	23.14289	24.7352	VRAI	6.41806	32.35711	123.14331	156.58185	475.78322	552.86328	64.56822	-0.50195
44.01	16	VRAI	6.50302	33.21388	123.57576	167.88421	521.63817	554.20428	23.292	24.71895	VRAI	6.4342	32.41862	123.3398	156.6303	475.86959	552.0385	64.6434	-0.50195
46.01	17	VRAI	6.51561	33.34262	123.78977	168.32939	521.07953	554.43903	23.80062	24.7392	VRAI	6.4051	32.22123	129.1232	155.95444	475.36959	551.64697	64.85546	-0.50195
48	17	VRAI	6.51891	33.356	123.73235	168.48146	521.47028	554.0708	23.22879	24.72374	VRAI	6.4215	32.38842	128.78401	157.1741	475.85559	552.1672	64.4821	-0.50195
50.02	17	VRAI	6.5191	33.28262	123.45886	168.58743	521.33891	553.87705	23.31759	24.71894	VRAI	6.42766	32.44554	128.14255	156.67442	476.02503	551.87328	64.58748	-0.50195
52	18	VRAI	6.5287	33.31254	123.78258	168.18727	521.4959	554.28662	23.3545	24.7392	VRAI	6.42173	32.38871	128.71862	157.2716	475.75883	552.63084	64.86676	-0.50195
54	19	VRAI	6.50228	33.20532	123.58033	168.65969	521.69136	553.86737	22.9573	24.71895	VRAI	6.43383	32.51686	128.52125	156.75523	476.23778	551.65056	64.5717	-0.50195
56.02	19	VRAI	6.50417	33.21383	123.37398	168.58039	521.15443	554.07666	23.7089	24.70817	VRAI	6.42186	32.38009	128.54484	156.26554	475.7238	551.14451	64.54768	-0.50195
58	20	VRAI	6.52286	33.1493	124.32182	167.95934	521.52557	554.24581	23.3371	24.72763	VRAI	6.42585	32.43027	123.27265	156.71988	475.97232	552.12555	64.77223	-0.50195
60	21	VRAI	6.505	33.23414	123.71851	168.63088	521.93714	553.82357	22.89786	24.75541	VRAI	6.41012	32.30394	128.51034	157.04828	475.62884	552.25236	64.94459	-0.50195
62.02	22	VRAI	6.51046	33.32083	123.7775	168.2489	521.51835	554.1001	23.2636	24.7393	VRAI	6.43091	32.4814	129.23432	157.01016	475.74775	551.94855	64.53756	-0.50195
64.02	23	VRAI	6.51703	33.35717	124.05712	168.10959	521.34188	554.430	23.62498	24.73541	VRAI	6.42267	32.39822	128.23006	156.68889	475.92734	552.30536	64.48079	-0.50195
66	23	VRAI	6.50277	33.3073	123.64548	167.8974	521.54883	554.11053	23.23973	24.70817	VRAI	6.43383	32.43071	123.18005	156.68441	476.3853	551.82808	64.12945	-0.50195
68	23	VRAI	6.53188	33.50737	124.70319	167.52751	520.86407	553.8894	23.5237	24.7392	VRAI	6.40664	32.23667	128.82828	156.63676	475.5488	551.76575	64.72283	-0.50195
70.01	23	VRAI	6.50176	33.2002	123.46574	168.06788	521.53551	554.10736	23.2332	24.71894	VRAI	6.42721	32.44488	128.06281	157.11603	476.01253	552.42326	64.4038	-0.50195
72	24	VRAI	6.53076	33.43978	124.12445	168.66573	521.12328	554.24218	23.58986	24.75541	VRAI	6.40861	32.29957	123.55884	156.88819	475.81223	551.88914	64.48138	-0.50195
74.01	25	VRAI	6.50359	33.28171	123.83142	167.33721	521.82085	553.98309	22.85884	24.72374	VRAI	6.40579	32.2281	128.51048	156.73889	475.62239	552.05891	64.70783	-0.50195
76	27	VRAI	6.53066	33.43669	124.1734	168.339	520.89978	553.7724	23.51777	24.72374	VRAI	6.41024	32.30311	128.33995	156.60449	475.96783	551.68671	64.2824	-0.50195

6 CLINICAL APPLICATION

6.1 Training requirements

The users should be familiar with Ophthalmology. It is important to provide the patients with the clinical judgment from the side of a Physician, to verify the suggested protocols according to the specific patient under consideration for treatment.



CAUTION: Device shall be operated only by medical staff with experience in Ophthalmology.

According to their expertise, such professionals can decide the correct use of the Device depending on the type of application

6.2 Intended use

RetinaWise measures the patient's pupil diameter as a function of time and specific scotopic light stimuli of the retina.

Disclaimer:

Eye pupil diameter variation depends on the level of light exposure. Eye pupil dynamics and pupil reflex is a patient's eye health index. The physician is the sole responsible to define any diagnosis or therapy based on the measurement taken with *retinaWISE*.

6.3 Patient Population

No restriction.

6.4 Contraindications

Patient with pathology as Migraine, Uveitis, Albinism, Aniridia Blepharitis, Coloboma, Amblyopia, and other need to be evaluated by medical doctors before being subjected to the analyses.

6.5 Side Effects

No side effect.

6.6 Residual Risk

No residual risks.

6.7 Warnings related to the applications

No warnings related.

Contact your Customer Service Representative for further instructions.

6.8 Precautions

Precaution shall be taken during Device operation since it requires dark or extremely dimmed light examination rooms. In such conditions patients might be disoriented. The room's floor must be free of obstacles to avoid accidents.

The examination could take several minutes up to 10 minutes. It is recommended to configure the exam room to be quiet, with a pleasant controlled ambient temperature so as to allow the patient to relax for optimal data acquisition.

6.9 Adverse events management

In case of adverse events such as.

- fuse failure or repeated fuse failure,
- stimulation spectroscopy does not correspond to specifications,
- field of view is not homogeneous,

manage them according to appropriate clinical practices and inform us (see §9.3 for contacts).

7 MAINTENANCE

This Device is designed for maximum safety and performance.

Ordinary & Extraordinary Maintenance of the Device	
Under careful and normal operating condition	Intensive use, dust, or continuous relocation
Every 6 Months	More frequent review required: every 3 months

7.1 Checking the line cable

The Device is provided with an external power supply. The line cable of the power supply can deteriorate over time and therefore it is necessary to periodically check the status of the line cable. In case of damages, contact Oculox Technologies for a replacement power supply.



Warning: Power supply cord can be replaced by authorized personnel only

7.2 Checking the labels

It is under the user's responsibility to keep the safety labels in good condition.
It is necessary to replace all the labels that are damaged.

7.3 Cleaning

All the cleaning procedures must be performed after the Device is switched off.



Warning: Before performing any cleaning action, be certain that the Device is switched OFF

7.3.1 Cleaning of the "Headrest"

Applied parts such as the forehead band, the chin rest cup and the handles of the headrest are made of easy-to-clean materials.

To comply with general hygiene requirements and to prevent the transmission of infections, the applied parts should be disinfected prior to every examination.

Several disinfection solutions exist, e.g., with 70% isopropyl alcohol, or ready-for-use disposable 70% ethanol disinfectant wipes.


Surface friendly disinfectants (containing aldehyde or aldehyde-free) are also permitted, such as Kohrs Olin FF.

**WARNING!**

- The efficacy of the disinfectants are determined by the disinfection manufacturers.
- The efficacy of the disinfectants mentioned above was not tested on their correct disinfection effect on the headrest.
- The efficacy of the disinfectants must be guaranteed by the user or the reprocessing responsible person with the validation of their own disinfection process.
- Comply with the stipulated exposure time.
- Observe the manufacturer's safety instructions.
- Too strong or aggressive disinfectants or cleaning liquids e.g. hydrogen peroxide will damage the finish and coating of the device.
- Do not use sprays.
- Do not use any cloths that drip.
- Wring out any soaked cloths before use when necessary.
- Ensure that no liquid penetrates the device.

7.3.2 Cleaning of the Cameras

Before cleaning the lenses, verify that there are not dust on the surface which can damage the lenses surface, remove the dust with a Optic Bulb Blowers or a dust and oil free air flow.



Once removed the dusts, the surface can be cleaned with Isopropanol and a Microfiber or equivalent tissue.

7.3.3 Cleaning of the external covering panels

Clean periodically the external covering panels of the Device using a cloth dampened with a standard cleaning solution. Avoid spraying directly the cleaning solution on the panels.



Warning: Avoid penetration of the cleaning solution inside the Device.

7.3.4 Contact for Information Request

If any further information is needed about cleaning, please contact Oculox Technologies

8 TROUBLESHOOTING

Issue	Possible cause	Solution
<i>The device does not switch ON</i>	The power cord is not correctly inserted in the grid wall plug	Control cabling
<i>Device does not connect to PC</i>	Device is not recognised by PC's USB port	Restart PC to reset and allow USB port configuration
<i>Connection with monitoring camera fails</i>	System configuration is corrupted	Setup and configure system
<i>Monitoring camera image stream is not centred</i>	System configuration is corrupted	Setup and configure system
<i>Device joystick is blocked</i>	Joystick brake is tight and block the system	Unblock joystick brake
<i>Patient alignment fails</i>	Alignment procedure is not followed correctly	Refer to User Manual to properly align device

9 CUSTOMER SERVICE

9.1 Warranty and Manufacturer's responsibilities

The Limited Warranty provided by Oculox Technologies applies to physical goods, and only for physical goods (the "Physical Goods").

This Limited Warranty covers any defects in material or workmanship under normal use during the Warranty Period.

During the Warranty Period, Oculox Technologies will repair or replace, at no charge, products or parts of a product that proves defective because of improper material or workmanship, under normal use and maintenance. Oculox Technologies will either repair the Product at no charge, using new or refurbished replacement parts.

The Warranty Period for Physical Goods purchased from Oculox Technologies is **1 year (366 days) from the date of purchase**.

A replacement Physical Good or part assumes the remaining warranty of the original Physical Good or 366 days from the date of replacement or repair, whichever is longer.

This Limited Warranty does not cover any problem that is caused by:

- malfunctions or damage not resulting from defects in material or workmanship.

To obtain warranty service, you must first contact us to determine the problem and the most appropriate solution for you.

9.2 Repairs and modifications of the Device

- Only authorised service personnel can execute repairs and maintenance,
- It is recommended that you follow the standard maintenance program,
- It is recommended that you replace all damaged components,
- Use only original spare parts,
- Constructive modifications are not permitted.

9.3 Contacts

Oculox provides its customers with email and telephone troubleshooting, apart from the on-site customer trainings, repairs and maintenance.

Please contact Oculox Service Department through the contacts reported here below. Please keep on hand the Serial Number of your Device.

Technical Support
<ul style="list-style-type: none">• <u>Direct connection to technical support portal:</u> From the menu bar <i>Help</i> of RetinaWISE programm, selection of <i>Technical support</i>, and this will open a direct web link• <u>From the link:</u> https://oculox.atlassian.net/servicedesk/customer/portal/3

9.4 Reporting of serious incidents

Any serious incident that has occurred in relation to the Device should be reported to the Manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Service Department
<u>Oculox Technologies SA</u> Via Industria 3 CH-6933 Muzzano Switzerland @-mail: service@oculox.com Phone: +4191 210 89 61

10 TECHNICAL SPECIFICATIONS

10.1 General Specifications














Classification according to Medical Device Regulation	Class I	
LED Classification according to IEC 62471	Exempt	
Mains	Only supplied power supply can be used – Mean Well GSM90B24-P1M	
	Rating:	100-240 VAC – 50/60 Hz 24V DC
Max power absorption	40VA	
Type of protection against electric shock	Class II	
Applied part	Type B	
Degree of protection against electric shock (EDS)	Level 4	
Degree of protection against the ingress of liquid	IP X0	
Mode of operation	Continuous	
Fuse	5x20mm – 4AF	
Dimension	350 (D) x 330 (W) x 400 (H) mm	
Weight	10 kg	
Operating temperature	10°- 30° C	
Storage/Transport Temperature	Min +5° C / max +40° C	
Operating humidity	30% - 85% (non-condensing)	
Operating pressure	From 800 hPa to 1060 hPa	
Max noise (dBA)	<70 dB	
Shelf life	5 years	

10.2 LED specifications

LED type	General purpose LED
Wavelength	420, 450, 470, 520, 590 and 630 nm
Max irradiance	5 (W/m ²)

11 GLOSSARY

A Glossary of symbols and abbreviations used in this Manual is reported below.

Symbols	Description
	Read the enclosed documentation
	CE Mark
	Symbol of applied part type B, according to standard 60601-1
	Symbol indicating that the Device cannot be disposed of as municipal waste, but must be separated in accordance with the WEEE (Waste Electrical and Electronic Equipment)
	Manufacturing date
	Manufacturer
	Serial Number
	Medical Device Symbol
	Catalogue number
	Authorized representative in the EU
	UDI label
	Attention - consult internal documentation
	Power ON/OFF

12 Appendix A: EMC Tables



CAUTION!

To guarantee the safety of the user, the patient and others, use only accessories and spare parts specified by the Manufacturer of this product. Other accessories or spare parts can cause the emission of increased electromagnetic radiation or reduced immunity against interference.

IMPORTANT!

Medical electrical devices are subject to special precautions with regard to electromagnetic compatibility (EMC) according to IEC 60601-1-2. Make sure you observe the notes on EMC for installation and operation. Medical electrical devices can be influenced by mobile HF communication devices (i.e. mobile phone). If it is necessary to stack the devices or place them next to each other, and HF interference is observed, make sure you observe the intended use of the devices.

IEC 60601-1-2 Table 201		
Guidance and Manufacturer's declaration – electromagnetic emission		
RetinaWISE is intended for use in the electromagnetic environment specified below. The customer or the end user of RetinaWISE should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment
RF emission – CISPR 11	Group 1	RetinaWISE uses RF energy only for its internal function. As a result, its RF emissions are very low and are not likely to cause any interference in nearby electronic devices.
RF emission – CISPR 11	Class B	RetinaWISE is suitable for use in all environments including households and those directly connected public low-voltage grid supplying buildings used for domestic purposes.
Harmonic emission IEC 61000-3-2	Class A	
Voltage fluctuation/flicker emission IEC 61000-3-3	Compliant	

IEC 60601-1-2 Table 202			
Guidance and Manufacturer's declaration – electromagnetic immunity			
RetinaWISE is intended for use in the electromagnetic environment specified below. The customer or the end user of RetinaWISE should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment
Electrostatic discharge (ESD) IEC 61000-4-2	In contact ± 8 kV	In contact ± 8 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
	In air ± 2; 4; 8; 15 kV	In air ± 2; 4; 8; 15 kV	
Electrical Fast Transient/Burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
	±1 kV for I/O lines	±1 kV for I/O lines	
Surges IEC 61000-4-5	±1 kV differential mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
	±2 kV common mode	±2 kV common mode	
Voltage Dips, Short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U _T for 0.5 cycle	0% U _T for 0.5 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of RetinaWISE requires continuous operation during power mains interruptions, it is recommended that RetinaWISE be powered from an Uninterruptible Power Supply or Battery.
	0% U _T for 1 cycle	0% U _T for 1 cycle	
	70% U _T for 25 cycles	70% U _T for 25 cycles	
	0% U _T for 250 cycles	0% U _T for 250 cycles	
Magnetic field at mains frequency (50/60Hz) IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Proximity Magnetic Field IEC 61000-4-39	30 KHz, CW, 8 A/m	30 KHz, CW, 8 A/m	Magnetic fields generated by RFID devices should have levels characteristic of a typical location in a commercial or hospital environment.
	134.2 KHz, PM, 65 A/m	134.2 KHz, PM, 65 A/m	
	13560 KHz, PM, 7.5 A/m	13560 KHz, PM, 7.5 A/m	
Note: U _T is the AC mains voltage prior to application of the test level.			

IEC 60601-1-2 Table 204**Guidance and Manufacturer's declaration – electromagnetic immunity**

RetinaWISE is intended for use in the electromagnetic environment specified below. The customer or the end user of RetinaWISE should ensure that it is used in such an environment. Portable and mobile RF communication equipment should not be used within 30 cm of any part of the device including cables.

Immunity test	IEC 60601 Test level	Compliance level		Electromagnetic environment
Conducted RF IEC 61000-4-6	3 V _{eff} from 150 kHz to 80 MHz 6 V _{eff} in ISM frequency	3 V _{eff}		d = 30 cm
Radiated RF IEC 61000-4-3	3 V/m from 80 MHz to 2.7 GHz	3 V/m		d = 30 cm
Immunity to proximity fields from wireless RF communication devices	TETRA 400 380 – 390 MHz	27 V/m	27 V/m	d = 30 cm
	GMRS 460 FRS 460 430 – 170 MHz	28 V/m	28 V/m	
	LTE Band 13, 17 704 – 787 MHz	9 V/m	9 V/m	
	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5 800 960 MHz	28 V/m	28 V/m	
	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 5 1700 – 1990 MHz	28 V/m	28 V/m	
	Bluetooth, WLAN, 802.11 b/g/n, RIFD 2450, LTE Band 70 2400 – 2570 MHz	28 V/m	28 V/m	
	WLAN 802.11 a/n 5100 – 5800 MHz	9 V/m	9 V/m	

Note 1: at 80 MHz and 800 MHz, the higher frequency range applies

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Device is used exceeds the applicable RF compliance level above, the Device should be observed to verify normal operation.

If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Device.

b) Over the frequency range 150 KHz to 80 MHz, field strength should be less than 3 V/m

IEC 60601-1-2 Table 206**Recommended separation distances between portable and mobile RF communication equipment and the Device**

The equipment is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications

equipment (transmitters) and the Device as recommended below, according to the maximum power of communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 KHz to 80 MHz $d=1,17*\sqrt{P}$ m	80 MHz to 800 MHz $d=1,17*\sqrt{P}$ m	800 MHz to 2,5 GHz $d=2,33*\sqrt{P}$ m
0,01	0,117	0,117	0,233
0,1	0,370	0,370	0,740
1	1,17	1,17	2,33
10	3,70	3,70	7,40
100	11,7	11,7	23,3

For transmitters rated at maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer

Note 1: at 80 MHz and 800 MHz, the higher frequency range applies

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.



BURST:

Upon the occurrence of a voltage transient (BURST) in the device power supply, the communication with the cameras could be lost with a consequent error message. At the end of the disturbance, the SW might prompt again to load the program, then from the main switch turn off and on again.



Voltage dips/Voltage interruption:

In case of a voltage gap during normal operation, when power returns, the machine may restart by emitting continuous flashes. When power returns, to resume normal operation, permanently turn off the machine from the main ON/OFF switch and then restart.



ESD:

Upon the occurrence of a possible electrostatic discharge induced toward the machine, communication with the cameras could be lost, with a consequent error message. To resume normal operation, permanently turn off the machine from the main ON/OFF switch and then restart.